

mediven elegance®, mediven® for men, mediven active®, mediven® comfort, mediven plus®

Circumferences

Sizes	I	II	III	IV	V	VI	VII
eG thigh with silicone topband	43 – 48	45 – 52	49 – 56	53 – 60	56 – 64	60 – 68	64 – 72
eG thigh with silicone topband extra wide	49 – 57	53 – 62	57 – 67	61 – 72	65 – 77	69 – 81	73 – 85
eG	43 – 57	45 – 62	49 – 67	53 – 72	56 – 77	60 – 81	64 – 85
eF	39 – 52	41 – 56	44 – 60	47 – 65	50 – 69	53 – 73	57 – 78
eE	30 – 37	33 – 40	35 – 43	37 – 45	39 – 48	41 – 51	42 – 53
eD	27 – 33	29 – 36	32 – 39	34 – 42	36 – 45	38 – 48	40 – 50
eC	28 – 34	30 – 37	33 – 40	35 – 43	37 – 46	39 – 49	41 – 51
eB ₁	23 – 27	24 – 29	26 – 32	29 – 35	31 – 37	33 – 39	35 – 41
eB	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
eY	26 – 31	28 – 33	29 – 35	31 – 37	32 – 38	33 – 40	34 – 42
eA	17 – 22	19 – 24	21 – 26	23 – 29	25 – 32	27 – 34	28 – 35

Length

eG regular	72 – 83
eG petite	62 – 71
eD regular	39 – 44
eD petite	34 – 38

Length mediven for men

eD regular	39 – 44
eD petite	34 – 38


Length mediven active

eD regular	43 – 49
eD petite	37 – 42

mediven® forte

Circumferences

Sizes	I	II	III	IV	V	VI	VII
eG thigh with silicone topband	44 – 47	46 – 50	50 – 54	53 – 58	56 – 61	60 – 65	64 – 70
eG thigh with silicone topband extra wide	48 – 53	51 – 57	55 – 61	59 – 65	62 – 70	66 – 74	71 – 80
eG	44 – 53	46 – 57	50 – 61	53 – 65	56 – 70	60 – 74	64 – 80
eF	40 – 47	42 – 51	45 – 55	48 – 58	51 – 62	54 – 67	58 – 72
eE	30 – 37	33 – 40	36 – 42	38 – 44	40 – 46	42 – 48	44 – 50
eD	28 – 33	30 – 35	32 – 37	34 – 39	36 – 41	38 – 43	40 – 45
eC	29 – 34	31 – 36	33 – 38	35 – 40	37 – 42	39 – 44	41 – 46
eB ₁	24 – 26	26 – 28	28 – 31	30 – 33	32 – 35	34 – 37	36 – 39
eB	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
eY	28 – 30	29 – 32	31 – 34	33 – 36	34 – 38	36 – 40	37 – 41
eA	18 – 21	20 – 23	22 – 25	24 – 27	26 – 29	28 – 31	30 – 33



Length

eG regular	72 – 83
eG petite	62 – 71
eD regular	38 – 43
eD petite	33 – 37

Toe

Standard with open toe
Short foot sizes 36 – 38
Normal foot sizes 38 – 43