

The order cannot be processed without details! Option field	
	QUOTE
Invoice address (stamp or complete legibly) Customer No.	Shipping address (if different from invoice address)
Company	Company
Street	Street
ZIP/Place	ZIP/Place
Type of order	
Mounting by customer ²⁾ (Please fill out the "Wheelchair details" section completely!)	1
Mounting by Alber in the factory ³⁾ 1565790	Consignment
 twion T24 Scope of delivery: 2 drive wheels, 6 km/h, 24" wheel size Standard model: 6.25 kg per wheel 2 built-in battery packs in the wheels 1 charger 	
*Please refer to your price list	
Wheelchair details	
Wheelchair manufacturer (e.g. Küschall)	<u> </u>
Model (e.g. Compact)	
Year of construction	* Seat width (cost tube measured on the outside)
Variants (e.g. Hemi)	••••••••••••
Drum brakes ⁵⁾ Yes No Note: Drum brakes must be removed!	1
Equipped with the vehicle retention system (please consult Alber customer service)	Seat height (measured to lower edge of seat sling)
Equipped with rear wheel extension	Wheel combar (if people)
Available wheel position	Wheel camber (if possible)
1 = rear	□ 3°
2 = middle 123	⊔ °
3 = front	
Important accessories	-
Bracket to fit on the wheelchair	
(in a few cases price could be different)	
Anti-tippers with jack-up function 1489214 (Pair)	
Comments (mounting, delivery, etc.)	
	DIN EN ISO 9001

Ser√ice tested







Option field



The order cannot be processed without details!



Со			
	ver		
0	White cover		(Pair)
	Black cover		(Pair)
Du	sh-rim		
ru	SII-riiii		
0	Matt aluminium, coated black		(Pair)
	Stainless steel, brushed (+400 g per wheel)		(Pair)
	Carbolife CURVE L		(Pair)
Ту	res		
0	Schwalbe Marathon Plus Evolution		(Pair)
	Schwalbe One (-300 g per wheel)		(Pair)
Ac	cessories		
	Mobility Plus package For smartphones with the Android or iOS operating system and able to connect to the Comprising: Outer ring with exclusive design, b operating instructions, licence code Functions: Reading off driving parameters including batte speed, number of kilometres covered 2 additional ECO and SPORT modes Increase assist up to 10 km/h BlueDrive remote control using a smartphone Records journeys and trip profiles	ooklet with ery status,	
	Brake extension		
	Ergonomic stub axle (set)	1591500	
	Stub axle for transport	1591362	
	anual wheel axle retrofit 4) 5) (con Retrofit, stub axle to ø ½"	npatible with tw	vion wheels) (Pair)

- In order to ensure that your order is processed smoothly and promptly, we require ALL the information in the shaded fields.
 Mounting kits can only be sent to trained specialist dealers (mounting training).
 Please send the wheelchair carriage with a copy of this order as soon as possible and carriage-paid to Alber.
- If the manual wheels are used in combination with twion, these must be fitted with ½" stub axles, without drum brakes and 1" wide tyres. If drum brakes are installed, these must be removed as they are not compatible with twion drive wheels. Standard lever brakes are used instead.