

☒ Standard

☐ Option field

☐ **Required fields:**<sup>1</sup>

Information required for processing



☐ ORDER

☐ QUOTE

**Invoice address** (stamp or fill in legibly)

Kd.-Nr. \_\_\_\_\_

Company

Street

ZIP/Place

**Shipping address** (if different from invoice address)

Company

Street

ZIP/Place

### Type of order

☐ Self-fitting<sup>2</sup> (fully complete the section "Information about the wheelchair"!)

☐ Mounting at the Alber factory<sup>3</sup> 1592525

### Delivery

☒ Standard

0,00 EUR

### Consignment

#### ☐ e-motion DuoDrive

With control unit (bracket included) 1594352  
Wheel size 22" (550 mm)

#### ☐ e-motion DuoDrive

With control unit (bracket included) 1594353  
Wheel size 24" (600 mm)

#### ☐ e-motion DuoDrive

With control unit (bracket included) 1594355  
Wheel size 25"  
625 with Right Run tyres, 615 with Marathon Plus tyres

#### ☐ Upcharge ECS Remote Control 1592486

Rollback delay can only be activated with  
ECS Remote Control or smartphone app with Mobility Plus Package

☐ Extended guarantee Li-Ion battery pack 1592523  
(to 36 months)

### Important accessories

☐ Bracket for mounting to the wheelchair

☐ Anti-tippers with jack-up function 1593206 (set of two)

### Wheelchair details

Wheelchair manufacturer  
(e.g. Meyra)

Model (e.g. Eurochair, 1.850)

Year of manufacture

Version (e.g. Hemi)

Drum brakes<sup>5</sup>

☐ Yes

☐ No

↳ Attention: Drum brakes must be removed

☐ Car restraint system available (please consult with Alber customer service)

☐ Wheel base extension available

Existing wheel position

☐ 1 = back

☐ 2 = centre

☐ 3 = front



Seatwidth


(seat tube measured externally)



\_\_\_\_\_ cm

Seat height back

(measured to lower edge  
of seat sling)



\_\_\_\_\_ cm

Diameter footrest  
for bracket of control  
unit DuoDrive



Ø \_\_\_\_\_ mm

Wheel camber (if possible)

☐ 0°

☐ 3°

**Comments** (mounting, delivery etc.)

☒ Standard

☐ Option field

☐ **Required fields:**<sup>1</sup>

Information required for processing

### Tyres

☒ Schwalbe Right Run (Standard) black 1592530  
(1" for wheel size 22" and 24", 1,1" for wheel size 25")

☐ Schwalbe Marathon Plus Evolution 1565894  
with SmartGuard puncture protection, black  
(not available for 22")

☐ Pro Speed puncture-proof solid PU 1" 1592532  
made of solid PU, black  
(+ 300 g per wheel, not available for 22" and 25")

☐ Puncture-proof with inlay 1 3/8" 1592533  
with PU insert and tyre cover, schwarz  
(+ 500 g per wheel, not available for 25")

### Push rims

☒ Stainless steel (standard) 1592526

☐ Coated push rim Tetra 1592527  
(not available for 25")  
(+250 g per wheel)

☐ Carbolife CURVE L 1592528  
powder coated, black  
(not available for 22")

☐ Carbolife QUADRO 1592529  
Tetra-Gripp coating  
(Only available for 24")

### Push rim covers (retrofit)

☐ push rim cover 22" 1592425 (pair)

☐ push rim cover 24" 1592421 (pair)

☐ push rim cover 25" 1592426 (pair)

### Spoke cover (retrofit)

☐ Spoke cover 22", Design „Comic“ 1592266 (pair)

☐ Spoke cover 24", Design „Sphere“ 1592268 (pair)

☐ Spoke cover 25", Design „Sphere“ 1592269 (pair)

☐ Spoke cover 24", plain 1592553 1592553 (pair)

### Accessories

☐ Ergonomic stub axle 1592972 (pair)

☐ Transport Axle 1591362

☐ Conversion of manual wheels to 1/2" stub axle<sup>4</sup> 1565005

☐ Brake extension kit 1467780 (pair)

Date  /  /  Signature of ordering party

Date 01.05.2021. All prices before VAT, including packaging ex factory. The general Terms and Conditions of Alber GmbH apply. Subject to changes.

<sup>1</sup> In order to smoothly and quickly process your order we require ALL information in the fields highlighted in blue.  
<sup>2</sup> Installation kits can only be shipped to trained specialist dealers (installation training).  
<sup>3</sup> Please send the wheelchair with a copy of this order as soon as possible and freight paid to Alber.  
<sup>4</sup> If the manual wheels are used in combination with e-fix, they must be equipped with 1/2" stub axles, without drum brakes and tyres of size 1".  
<sup>5</sup> If present, drum brakes must be dismantled as they are not compatible with e-fix drive wheels. The use of standard lever brakes is required.