

Hepro Work Chairs



Table of Contents

- 3 Functionality
- 4 Product naming key
- 5 Hepro Tilto
- 8 Hepro Aktiv
- 10 Hepro Luna
- 12 Accessories overview
- 14 Accessory compatibility guide



Hepro Tilto



Hepro Aktiv



Hepro Luna

Functionality

Hepro work chairs are robust and functional. They are top quality and hard wearing. The chairs are safe, have readily accessible levers and feature many adjustment options for personalised sitting comfort.

This is an overview of the functions that feature in all our models.

Ergonomic backrest

Flexible Armrests

Regulation options

Large foot space

Robust wheels



Product naming key

In the various chair title's, we use different indicators to separate them. Here we are going to explain the different parts of the name with the Example of a chair out of the Tilto Series. Hepro G2 Tilto 15 EHR

EHR / HR – Size

The Tilto chairs have either HR or EHR, HR stands for Høy Rygg, which translates into High Backrest EHR stands for Ekstra Høy Rygg, which translates into Extra High Backrest

15 / 10 – Height

15 represents the height adjustment of the gas lift/ electric lift. In this case the gas lifted chair has a 15cm height adjustment. Alternatively 10 would mean a 10cm height adjustment.



Elektrik mechanism

G2 / E2 – Lift

Here G in G2 represents Gas lifted chair, and 2 means the chair has rotation/lockable rotation. E would represent Electric lift. And 1 would mean the chair has locked rotation



Hepro Tilto

Hepro Tilto has both negative and positive seat tilt and is suitable for those who wish to vary their sitting position frequently. The seat tilt is balanced according to body weight, with an adjustment knob under the seat. With the tilt lever in the open position, the chair has a rocking function. It can be locked steplessly in the desired position, and the individual back angle adjustment is standard.

Head- and backrest

A headrest and inflatable lumbar support are standard. Together with an extra high backrest, this provides a good resting position. All chairs are delivered with foldable armrests.

Tilt function

The unique feature of our Tilto models is the tilt function which contributes to variation and movement. When working, it is easy to end up in a forward-leaning, static sitting position that results in impaired circulation.

With the tilt function, in the open position, the chair will:

- follow the movements of the body - all the way back to the relaxed sitting position
- opening the airways and stretching out the upper body and aiding good blood circulation

Movement adds to variation, relieving pressure on vulnerable pressure points and aiding good blood circulation to the pelvis, hips and thighs.



Tilto 24/7

Tilto 24/7 is the chair for you if you sit in your chair for many hours a day over a longer period and have a greater need for relief and comfort. Extra consideration has therefore been given to ergonomics and sitting comfort.

The chair has extra cushioning in the seat and back-rest and can be upholstered in special 24/7 breathable material so that you can sit comfortably for longer periods. The model has many individual adjustment options to suit your needs. Like the other Tilto models, the chair is designed to activate inactive users by means of the tilt function. This model has inflatable lumbar support. This model has inflatable lumbar support.

Hepro Tilto models



Hepro E2 Tilto 10/15 HR
Article nr. 30072/30356

Hepro G2 Tilto 10/15 HR
Article nr. 30321/30323



Hepro E2 Tilto 10/15 EHR
Article nr. 30366/30358

Hepro G2 Tilto 10/15 EHR
Article nr. 30320/30322



Hepro E2 Tilto 10/15 EHR 24-7
Article nr. 30362/30360

Hepro G2 Tilto 10/15 HR 24-7
Article nr. 30363/30364

Specifications

| Standard configuration | E2 Tilto 10 HR | E2 Tilto 15 HR | G2 Tilto 10 HR | G2 Tilto 15 HR |
|------------------------|----------------|----------------|----------------|----------------|
| Seat width | 47 cm | 47 cm | 47 cm | 47 cm |
| Seat depth | 48 cm | 48 cm | 48 cm | 48 cm |
| Seat height | 52 - 62 cm | 57 - 72 cm | 44,5 - 54,5 cm | 49 - 63 cm |
| Backrest width | 44 cm | 44 cm | 44 cm | 44 cm |
| Backrest height | 51 cm | 51 cm | 51 cm | 51 cm |
| Height adjustment | 10 cm | 15 cm | 10 cm | 15 cm |
| Seat angle | -6° +10° | -6° +10° | -20° +5° | -20° +5° |
| Armrest length | 35 cm | 35 cm | 35 cm | 35 cm |
| Width between armrests | 36 - 52 cm | 36 - 52 cm | 36 - 52 cm | 36 - 52 cm |
| Chair maximum width | 68 cm | 68 cm | 68 cm | 68 cm |
| Max. user weight | 150 kg | 150 kg | 150 kg | 150 kg |
| Article no. | 30072 | 30356 | 30321 | 30323 |

| Standard configuration | E2 Tilto 10 EHR | E2 Tilto 15 EHR | G2 Tilto 10 EHR | G2 Tilto 15 EHR |
|------------------------|-----------------|-----------------|-----------------|-----------------|
| Seat width | 47 cm | 47 cm | 47 cm | 47 cm |
| Seat depth | 48 cm | 48 cm | 48 cm | 48 cm |
| Seat height | 52 - 62 cm | 57 - 72 cm | 44,5 - 54,5 cm | 49 - 63 cm |
| Backrest width | 44 cm | 44 cm | 44 cm | 44 cm |
| Backrest height | 62 cm | 62 cm | 62 cm | 62 cm |
| Height adjustment | 10 cm | 15 cm | 10 cm | 15 cm |
| Seat angle | -6° +10° | -6° +10° | -6° +10° | -6° +10° |
| Armrest length | 35 cm | 35 cm | 35 cm | 35 cm |
| Width between armrests | 36 - 52 cm | 36 - 52 cm | 36 - 52 cm | 36 - 52 cm |
| Chair maximum width | 68 cm | 68 cm | 68 cm | 68 cm |
| Max. user weight | 150 kg | 150 kg | 150 kg | 150 kg |
| Article no. | 30366 | 30358 | 30320 | 30322 |

| Standard configuration | E2 Tilto 10 EHR 24-7 | E2 Tilto 15 EHR 24-7 | G2 Tilto 10 EHR 24-7 | G2 Tilto 15 EHR 24-7 |
|------------------------|----------------------|----------------------|----------------------|----------------------|
| Seat width | 49 cm | 49 cm | 49 cm | 49 cm |
| Seat depth | 52 cm | 52 cm | 52 cm | 52 cm |
| Seat height | 52 - 62 cm | 57 - 72 cm | 44,5 - 54,5 cm | 49 - 63 cm |
| Backrest width | 44 cm | 44 cm | 44 cm | 44 cm |
| Backrest height | 62 cm | 62 cm | 62 cm | 62 cm |
| Height adjustment | 10 cm | 15 cm | 10 cm | 15 cm |
| Seat angle | -6° +10° | -6° +10° | -6° +10° | -6° +10° |
| Armrest length | 35 cm | 35 cm | 35 cm | 35 cm |
| Width between armrests | 36 - 52 cm | 36 - 52 cm | 36 - 52 cm | 36 - 52 cm |
| Chair maximum width | 68 cm | 68 cm | 68 cm | 68 cm |
| Max. user weight | 150 kg | 150 kg | 150 kg | 150 kg |
| Article no. | 30362 | 30360 | 30363 | 30364 |



Hepro G2 Aktiv is an ergonomic chair that is ideally suited to provide relief during everyday activities either at home or at the office.

The chair has been developed for optimal flexibility and mobility, with a gas stem lift and a rocking feature. Hepro G2 Aktiv has a plethora of individual adjustment options and features designed to accommodate various needs.

Hepro G2 Aktiv chair features:

- individual stepless back angle adjustment
- rocking feature
- floating tilt function
- rotation with locking in 4 positions
- foldable armrests
- available in medium and large sizes

The seat and back

The seat and backrest are ergonomically designed to provide improved seating comfort.

The seat and back can be angled simultaneously and adjusted steplessly. With the lever in the open position, the chair has a floating tilt. It can be locked steplessly in the desired position.

The rocking feature

Putting the rocking feature in the open position, allows the chair to track your movements to relieve pressure points and vary your sitting position. The rocking feature also makes it easier to rise up from the ergonomic chair.



Hepro Aktiv Models



Hepro G2 Aktiv 15 M
Medium
Article nr. 30312

Hepro G2 Aktiv 15 L
Large
Article nr. 30313

Specifications

| Standard configuration | G2 Aktiv 15 M | G2 Aktiv 15 L |
|------------------------|---------------|---------------|
| Seat width | 44 cm | 46 cm |
| Seat depth | 43 cm | 46 cm |
| Seat height | 47 cm | 47 cm |
| Backrest width | 37 cm | 37 cm |
| Backrest height | 29 cm | 43 cm |
| Height adjustment | + 15 cm | + 15 cm |
| Seat angle | -6° +9° | -6° +9° |
| Armrest length | 35 cm | 35 cm |
| Width between armrests | 43 - 52 cm | 43 - 52 cm |
| Chair maximum width | 68 cm | 68 cm |
| Max. user weight | 150 kg | 150 kg |
| Article no. | 30312 | 30313 |

Hepro Luna

Hepro Luna is a work chair designed to provide pressure relief, lower back support, and comfort, so that you can continue working in an office environment or at home.

- The back is designed to ensure support, but also to provide comfort.
- The chair incorporates excellent features and individual controls, taking care of the varied needs of the user, while providing flexibility and mobility.

Luna G2 (gas lift)

With Hepro Luna you get uncompromising ergonomics combined with innovation and design. Hepro Luna is a work chair with:

- manual (gas) seat lift
- tilt function
- ergonomic seat

Luna E2 (electric lift)

Luna E2 is an easy-to-handle work chair, ergonomically adjusted to the body. This model has an electric lifting and lowering function and is especially suitable for:

- users with poor balance
- decreased strength in legs or hands

Tilt Function

The tilt function contributes to variation of sitting positions and movement. In a work situation it is easy to end up in a forward and static sitting position leading to reduced blood circulation in the body. With the tilt function in the open position, the chair will follow the movements of the body, all the way back to the relaxed sitting position, opening the airways and stretching out the upper body, and aiding good blood circulation. Movement also contributes to variation in sitting position, relieving pressure on the exposed pressure points and helps good blood flow to the pelvis, hips and thighs.



Hepro Luna Models



Hepro G2 Luna Tilt Small
Small
Article nr. 32516

Hepro G2 Luna Tilt Large
Large
Article nr. 32517



Hepro E2 Luna Tilt Small
Small
Article nr. 32518

Hepro E2 Luna Tilt Large
Large
Article nr. 32519

Specifications

| Standard configuration | Hepro G2 Luna Tilt Small | Hepro G2 Luna Tilt Large | Hepro E2 Luna Tilt Small | Hepro E2 Luna Tilt Large |
|------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| Seat width | 45 cm | 45 cm | 45 cm | 45 cm |
| Seat depth | 42 cm | 42 cm | 42 cm | 42 cm |
| Seat height | 45 - 60 cm | 45 - 60 cm | 50 - 65 cm | 50 - 65 cm |
| Backrest width | 45 cm | 47 cm | 45 cm | 47 cm |
| Backrest height | 42 cm | 45 cm | 42 cm | 45 cm |
| Backrest angle | -2° +22° | -2° +22° | -2° +22° | -2° +22° |
| Seat angle | -7° +6° | -7° +6° | -7° +6° | -7° +6° |
| Armrest length | 35 cm | 35 cm | 35 cm | 35 cm |
| Width between armrests | 45 - 60 cm | 45 - 60 cm | 45 - 60 cm | 45 - 60 cm |
| Chair maximum width | 77 cm | 77 cm | 77 cm | 77 cm |
| Max. user weight | 165 kg | 165 kg | 165 kg | 165 kg |
| Article no. | 32516 | 32517 | 32518 | 32519 |

Accessories




| Product name | Art. nr. |
|-----------------------|----------|
| Armrest Comfort right | 31593 |
| Armrest Comfort left | 31594 |

Armrest Comfort is a slimmer and more stable armrest compared to its counterpart Armrest Royal (standard). The armrest pad can be flipped up to clear more space on the sides.



| Product name | Art. nr. |
|------------------------------|----------|
| Armrest Comfort right x-long | 31596 |
| Armrest Comfort left x-long | 31597 |

Armrest comfort x-long is similar to the regular comfort armrest in size, but the connection point is 6cm further back, which is more beneficial for use cases where a user needs to move sideways out of the chair, or to have easier access to the armrests in a more back tilted position.



| Product name | Art. nr. |
|-------------------|----------|
| Breaklever x-long | 32715 |

A 15 cm longer break lever.



| Product name | Art. nr. |
|-------------------|----------|
| Electrical brakes | 31674 |

The electrical break is an accessory for electric chairs, for activating the breaks using an electric switch instead of the regular handheld break.



Armrest
Comfort



| Product name | Art. nr. |
|-----------------------------------|----------|
| Basic/ Luna, Common footrest assy | 300-1147 |

An attachable footrest to the wheelbase. Which increases the size of the platform for your feet.

Accessory compatibility guide

| | G2 Aktiv 15 M | G2 Aktiv 15 L | G2 Luna Tilt S | G2 Luna Tilt L | E2 Luna Tilt S | E2 Luna Tilt L | E2 Tilto 10 HR | E2 Tilto 15 HR | G2 Tilto 10 HR | G2 Tilto 15 HR | E2 Tilto 10 EHR | E2 Tilto 15 EHR | G2 Tilto 10 EHR | G2 Tilto 15 EHR | E2 Tilto 10 EHR 24-7 | E2 Tilto 15 EHR 24-7 | G2 Tilto 10 EHR 24-7 | G2 Tilto 15 EHR 24-7 |
|--------------------------------------|------------------|------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|--------------------|--------------------|--------------------|--------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| Armrest Comfort right | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x |
| Armrest Comfort left | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x |
| Armrest Comfort right x-long | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x |
| Armrest Comfort left x-long | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x |
| Breaklever x-long | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x |
| Electrical brakes | | | | | x | x | x | x | | | x | x | | | x | x | | |
| Basic/ Luna, Common footrest assy | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x | x |

Spare part

A housing component to combine standard gas springs with Hepro wheelbase, without the need for lockable rotation.

| Product name | Art. nr. |
|--------------|----------|
| Cone pipe | 31545 |



Active and safe in all stages of life

Hepro was established in 1987 and has its headquarters and production facilities surrounded by majestic mountains in Rognan, by the Sjøfjord. We have grown to become a leading player in technical aids and welfare technology solutions, with offices spread across Norway. In close collaboration with occupational therapists, we develop products and solutions tailored to our customers' needs for support with mobility, balance, comfort, and hygiene. The goal is to provide safety, independence, and increased activity through functional and user-oriented products and solutions. Hepro began its journey with the development and production of aids for adults. Our expertise in this field has expanded, and today we offer a wide range of products. Our portfolio also includes Krabat, a brand dedicated to innovative aids for children. Krabat products

give children the opportunity to achieve mastery and progress in a fun and engaging way. At a time with a growing elderly population and limited healthcare resources, Hepro works to free up time and improve care through technological solutions. We take a holistic approach to welfare technology. Our own platform, Hepro Respons, scales extremely well and makes it possible to integrate new technologies and develop unique services for municipalities both in Norway and abroad. Hepro is one of Norway's leading companies in assistive devices and welfare technology, but our presence extends far beyond the country's borders. Our products and solutions are exported worldwide. Visit our website for contact information and find out how Hepro can help you create an active and safe everyday life.