

## Active wheelchair *Avantgarde DV.*

### Master order form

#### Customer

Cust. no.

Name

Street no.

PC, city

#### Shipping address (if different from customer address)

Name

Street no.

PC, city

Com.

Name \_\_\_\_\_ Date \_\_\_\_\_ Signature \_\_\_\_\_



**Cost estimate**



**Order**



**Attachment for custom fabrication enclosed**



By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

#### Standard equipment

Article number: 480F160=0\_AA03\_C

**Price: x.xxx,- €**

valid from 01.01.2024 - 31.12.2024

- Aluminium tube frame with single and double cross brace
- Choice of front frame offset
- Aluminium leg supports, individual leg supports, swing away, removable
- Adaptable back support upholstery
- Aluminium side panels
- Modular design
- Drive wheels with quick-release axles
- Knee lever wheel lock 4.0
- Calf strap

#### Folding, max. user weight 140 kg

Crash-tested according to DIN ISO 7176-19 (with single and double cross brace)

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.



# Active wheelchair *Avantgarde DV*.

## Master order form

Connection of power add-on drive		Price
<input type="radio"/> MZ23	Klaxon Twist pull device (only valued for Germany, Switzerland and France)	Add. charge
<input type="radio"/> MZ24	Klaxon Twist push device (only valued for Germany, Switzerland and France)	Add. charge
<input type="radio"/> MZ25	Klaxon Twist pull and push device (only valued for Germany, Switzerland and France)	Add. charge
<input type="radio"/> MZ26	Klaxon Klick	Add. charge

Front frame type		Price
<input checked="" type="checkbox"/> MA000	DV	Add. charge

Front frame offset		Price
<input type="radio"/> MA14	0 cm offset	–
<input type="radio"/> MA15	2 cm offset	–

Equipment		Price
<input type="radio"/> MA110	Standard (for using in wheelchair accessible vehicles according DIN ISO 7176-19)	–
<input type="radio"/> MA111	CLT-Version (incl.: MA114, MA100, MG15, MG200, MH05, ME04)	Add. charge
<input type="radio"/> MA112	CLT-Version, Ultra (incl.: MA114, MA100, MG15, MG200, MH05, MD18, ME04, MC18)	Add. charge

Design of drive wheel adapter		Price
<input type="radio"/> MA113	With different setting options (for using in wheelchair accessible vehicles according DIN ISO 7176-19)	–
<input type="radio"/> MA114	Permanently welded	Add. charge

Cross brace		Price
<input type="radio"/> MA100	Single, CLT (max. load capacity 100 kg, up to SW 42)	Add. charge
<input type="radio"/> MA101	Double	–

Type of wheelchair		Price
<input type="radio"/>	Standard	Add. charge
<input type="radio"/> MA50	Foot-propelled wheelchair (reduced ground clearance in delivery condition for hemiplegia fittings)*	–

\* Risk due to getting caught on obstacles. Please observe safety information in the instructions for use!

Seat width (SW) in cm										
<input type="radio"/> 32	<input type="radio"/> 34	<input type="radio"/> 36	<input type="radio"/> 38	<input type="radio"/> 40	<input type="radio"/> 42	<input type="radio"/> 44	<input type="radio"/> 46	<input type="radio"/> 48	<input type="radio"/> 50	<input type="radio"/> 52

Seat depth (SD) in cm									
<input type="radio"/> 36	<input type="radio"/> 38	<input type="radio"/> 40	<input type="radio"/> 42	<input type="radio"/> 44	<input type="radio"/> 46	<input type="radio"/> 48	<input type="radio"/> 50	<input type="radio"/> 52	<input type="radio"/> 54

Anterior seat height in cm								
<input type="radio"/> 38	<input type="radio"/> 39	<input type="radio"/> 40	<input type="radio"/> 41	<input type="radio"/> 42	<input type="radio"/> 43	<input type="radio"/> 44	<input type="radio"/> 45	<input type="radio"/> 46
<input type="radio"/> 47	<input type="radio"/> 48	<input type="radio"/> 49	<input type="radio"/> 50	<input type="radio"/> 51	<input type="radio"/> 52	<input type="radio"/> 53	<input type="radio"/> 54	

Frame accessories		Price
<input type="radio"/> MA50	Transport wheels, pair	Add. charge
<input type="radio"/> MA51	Tip-assist, right	Add. charge
<input type="radio"/> MA52	Tip-assist, left	Add. charge
<input type="radio"/> MA55	Anti-tipper, right	Add. charge
<input type="radio"/> MA56	Anti-tipper, left	Add. charge
<input type="radio"/> MA57	Crutch holder with tip-assist, right	Add. charge
<input type="radio"/> MA58	Crutch holder with tip-assist, left	Add. charge



# Active wheelchair *Avantgarde DV.*

## Master order form

Caster wheels		Price
<input type="radio"/> MF01	3" skater role	–
<input type="radio"/> MF30	4" solid rubber, black	–
<input type="radio"/> MF35	4" soft caster, aluminium	Add. charge
<input type="radio"/> MF31	5" solid rubber, black	–
<input type="radio"/> MF36	5" soft caster, aluminium	Add. charge
<input type="radio"/> MF04	5.5" soft caster, grey	–
<input type="radio"/> MF05	6" pneumatic, grey	–
<input type="radio"/> MF32	6" solid rubber, black	–
<input type="radio"/> MF37	6" soft caster, aluminium	Add. charge

Caster axle	
<input checked="" type="checkbox"/> MF90	Caster forks with threaded axles

Posterior seat height in cm							
<input type="radio"/> 36	<input type="radio"/> 37	<input type="radio"/> 38	<input type="radio"/> 39	<input type="radio"/> 40	<input type="radio"/> 41	<input type="radio"/> 42	<input type="radio"/> 43
<input type="radio"/> 44	<input type="radio"/> 45	<input type="radio"/> 46	<input type="radio"/> 47	<input type="radio"/> 48	<input type="radio"/> 49	<input type="radio"/> 50	<input type="radio"/> 51

Drive wheel size		Price
<input type="radio"/> MG01	22"	–
<input type="radio"/> MG02	24"	–
<input type="radio"/> MG00	25"	–
<input type="radio"/> SBK13335.000	26"	Add. charge

Drive wheel camber			
<input type="radio"/> 0°	<input type="radio"/> 2°	<input type="radio"/> 3°	<input type="radio"/> 4°




# Active wheelchair *Avantgarde DV.*

## Master order form


**Horizontal drive wheel position** (for frame version adjustable)

☐ A ☐ B ☐ C ☐ D ☐ E ☐ F ☐ G

From position E, the drive wheel adapter must be turned around  
– right side to the left side and left side to the right side.



**Drive wheel positions with CLT equipment**





# Active wheelchair *Avantgarde DV.*

## Master order form

Drive wheel rim		22"	24"	25"	Price
<input type="radio"/> MG06	Hollow rims	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-
<input type="radio"/> MG10	Spinergy SPOX		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG11	Spinergy LX		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG12	Spinergy LXL Blade: black spokes		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG12	Spinergy LXL Blade: white spokes		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG15	Infinity Ultralight (not possible in combination with MG 34)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG150	Infinity Ultralight, coloured (not possible in comb. with MG 34)		<input checked="" type="checkbox"/>		Add. charge
Colour set Infinity Ultralight/caster fork		22"	24"	25"	Price
<input type="radio"/>	Sparkle light orange		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/>	Marine blue		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/>	Candy red		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/>	Sparkle granny-smith		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG80	With drum brake / attendant*	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG83	One-arm drive with double handrim left		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG83	One-arm drive with double handrim right		<input checked="" type="checkbox"/>		Add. charge

\*: not possible in combination with camber

Tyres		22"	24"	25"	Price
<input type="radio"/> MG30	Pneumatic tyre 1 3/8"		<input checked="" type="checkbox"/>		-
<input type="radio"/> MG31	Schwalbe RightRun 1"	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-
<input type="radio"/> MG34	PU 1" (not possible in combination with MG 15/MG 150)		<input checked="" type="checkbox"/>		-
<input type="radio"/> MG341	PU 1" black (not possible in combination with MG 15/MG 150)		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG37	Schwalbe Marathon Plus Evolution 1"	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG371	Schwalbe Marathon Plus Evolution 1.1"		<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG200	Schwalbe ONE 1"		<input checked="" type="checkbox"/>		Add. charge

Handrims		22"	24"	25"	Price
<input type="radio"/> MG50	Aluminium, black coated	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-
<input type="radio"/> MG52	Aluminium, silver anodised, high	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG53	Aluminium, silver anodised, low	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG54	rubber-coated	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG56	Stainless steel, polished, high	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG57	Stainless steel, polished, low	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Add. charge
<input type="radio"/> MG63	Curve L, anodised		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG64	Curve L, Grip for users with limited gripping function		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge
<input type="radio"/> MG65	GEKKO, ergonomic handrim		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Add. charge

Handrim mounting		Price
<input type="radio"/> SBK13379.001	Mounting position more narrow to customer requirements <input type="text" value=""/>	Add. charge
<input type="radio"/> MG70	Mounting position narrow: 15 mm	-
<input type="radio"/> MG71	Mounting position wide: 25 mm	-

Wheelbase extension including wheel lock adapter (not possible in combination with shock absorber)		Price
<input type="radio"/> MG86	Extra long	Add. charge
<input type="radio"/> MG82	Extra long with second fitting	Add. charge
<input type="radio"/> MG46	Without	-



# Active wheelchair *Avantgarde DV.*

## Master order form

Quick-release axle system		Price
<input type="radio"/> MG27	Quick-release axle	-
<input type="radio"/> MG28	Quick-release axle with Tetra release	Add. charge

Drive wheel accessories		Price
<input type="checkbox"/> MG87	High-pressure air pump, 14 bar	Add. charge
<input type="checkbox"/> MG88	Airman-pump	Add. charge

Shock absorber		Price
<input type="radio"/> MG97	Damping "soft" (up to 60 kg)	-
<input type="radio"/> MG98	Damping "medium" (61 kg – 100 kg)	Add. charge
<input type="radio"/> MG99	Damping "hard" (101 kg – 140 kg)	Add. charge

Wheel locks		Price
<input type="radio"/> MH01	Knee lever wheel lock 4.0	-
<input type="radio"/> MH30	Knee lever wheel lock 4.0, long lever	Add. charge
<input type="radio"/> MH05	Scissor wheel lock	Add. charge
<input type="radio"/> MH06	Knee lever wheel lock for user and attendant	Add. charge
<input type="radio"/> MH08	Knee lever wheel lock for one-handed operation, right	Add. charge
<input type="radio"/> MH09	Knee lever wheel lock for one-handed operation, left	Add. charge

Wheel lock accessories		Price
<input type="radio"/> MH10	Wheel lock lever extension, fixed	Add. charge
<input type="radio"/> MH21	Wheel lock lever extension, pluggable	Add. charge

**Back support height in cm**

☐ 25    ☐ 27,5    ☐ 30    ☐ 32,5    ☐ 35    ☐ 37,5    ☐ 40    ☐ 42,5    ☐ 45    ☐ 47,5  
☐ 50    ☐ 52,5    ☐ 55

Back support system in cm		Price
<input type="radio"/> MD01	Back support, height-adjustable, straight <input type="radio"/> RH 25 - 32,5 cm <input type="radio"/> RH 32,5 - 40 cm <input type="radio"/> RH 42,5 - 50 cm <input type="radio"/> RH 52,5 - 55 cm	-
<input type="radio"/> MD09	Back support, height-adjustable, with lumbar curve <input type="radio"/> RH 35 - 40 cm <input type="radio"/> RH 40 - 45 cm <input type="radio"/> RH 45 - 50 cm <input type="radio"/> RH 50 - 55 cm	-
<input type="radio"/> MD10	Back support, angle-adjustable in 6° increments (-15° bis +9°) <input type="radio"/> RH 30 - 40 cm <input type="radio"/> RH 35 - 45 cm <input type="radio"/> RH 40 - 50 cm	Add. charge
<input type="radio"/> MD12	Firmly welded back support (not approved for use in wheelchair accessible vehicles) <input type="radio"/> 0° <input type="radio"/> -5° <input type="radio"/> -10° <input type="radio"/> -15°	Add. charge

Back support upholstery		Price
<input type="radio"/> MD03	Adaptable, padded	Add. charge
<input type="radio"/> MD05	Adaptable, breathable	-
<input type="radio"/> MD18	Ultra light	Add. charge
<input type="radio"/> MD06	Without	-

Back support pad widening		Price
<input type="radio"/> MD80	+ 2 cm	Add. charge
<input type="radio"/> MD81	+ 4 cm	Add. charge
<input type="radio"/> MD82	Without	-



# Active wheelchair *Avantgarde DV.*

## Master order form

Use of BAXX back support system (include order form, model-specific restrictions are possible, take into account the reduction of the seat depth by using the Baxx back support system)		Price
<input type="radio"/> MD44	Baxx aluminium back support (enclosed)	Add. charge

Push handles		Price
<input type="radio"/> MD50	Short	–
<input type="radio"/> MD51	Long	–
<input type="radio"/> MD52	Telescoping, long	Add. charge
<input type="radio"/> MD55	Height-adjustable/removable	Add. charge
<input type="radio"/> MD45	Height-adjustable/removable for Baxx	Add. charge
<input type="radio"/> MD56	Fold-down (max. load capacity 125 kg/275 lbs)	Add. charge
<input type="radio"/> MD59	Without	–

Back support accessories		Price
<input type="radio"/> MD58	Stabiliser bar	–

Side panels		Price
<input type="radio"/> ME01	With protection against the cold	–
<input type="radio"/> ME02	With splash guard and protection against the cold	Add. charge
<input type="radio"/> ME04	Carbon, with splash guard	Add. charge
<input type="radio"/> ME05	Desk side panel with arm support, short	Add. charge
<input type="radio"/> ME06	Desk side panel with arm support, long	Add. charge
<input type="radio"/> ME17	Side panel with arm support long, depth-adjustable	Add. charge
<input type="radio"/> SBK24811.010	Side panel with arm support, pluggable (not in combination with single cross brace CLT (MA 100))	Add. charge

Side panel accessories		Price
<input type="radio"/> ME50	Arm support padded, swing away (not possible in combin. with desk side panel)	–
	<input type="radio"/> short <input type="radio"/> long <input type="radio"/> long, extra high	
<input type="radio"/> ME60	Tray	Add. charge
<input type="radio"/> ME66	Tray, half side	Add. charge
	<input type="radio"/> left <input type="radio"/> right	
<input type="radio"/> ME61	Swivel unit for channel forearm pads, right *	Add. charge
<input type="radio"/> ME62	Swivel unit for channel forearm pads, left *	Add. charge

\*: replaces the existing arm support

Channel forearm pad		Price
<input type="radio"/> ME70	Channel forearm pad, one-piece, left	Add. charge
<input type="radio"/> ME71	Channel forearm pad, one-piece, right	Add. charge

Seat upholstery		Price
<input type="radio"/> MC15	Seat upholstery with storage compartment (up to SW 38 cm = 1 compartment; from SW 40 cm = 2 compartments)	Add. charge
<input type="radio"/> MC16	Upholstery readjustable	–
<input type="radio"/> MC18	Ultra light	Add. charge
<input type="radio"/> MC06	Without	–



# Active wheelchair *Avantgarde DV.*

## Master order form

### Seat cushions

		Price
<input type="radio"/> MC01	<b>Standard</b> (cushion thickness 5 cm)	Add. charge
<input type="radio"/> MC02	<b>Cubic Foam</b> (cushion thickness 5 cm)	Add. charge
<input type="radio"/> MC05	<b>Evolight</b> (cushion thickness 6 cm)	Add. charge
<input type="radio"/> MC20	<b>Terra</b> (cushion thickness 10 cm)	Add. charge
<input type="radio"/> MC21	<b>Terra Aquos</b> (cushion thickness 10 cm)	Add. charge
<input type="radio"/> MC22	<b>Terra Flair</b> (cushion thickness 10 cm)	Add. charge
<input type="radio"/> MC23	<b>Cloud</b> (cushion thickness 10 cm)	Add. charge
<input type="radio"/> MC29	<b>Livo</b> (cushion thickness 5 cm)	Add. charge
<input type="radio"/> MC30	<b>EvoFlair</b> (cushion thickness 5 cm)	Add. charge


Material	Aluminium (serrated, MB 160)	Plastic (MB 170)	Accessories
Leg support			
Single-panel, wide (MB 000)	<input type="radio"/> Standard (MB100) – <input type="radio"/> For short LL (MB150) AC	–	<input type="radio"/> With lock (MB70) AC –
Individual	<input type="radio"/> Standard 16 cm (MB010) –	<input type="radio"/> Standard (MB010) –	–
Elevating	–	<input type="radio"/> Pair (MB25) AC <input type="radio"/> Right* (MB26) AC <input type="radio"/> Left* (MB27) AC	–
Tube foot support	<input type="radio"/> Standard (MB120) –	–	<input type="radio"/> Plastic footplate (MB74) AC
<input type="radio"/> Without (MB80)	–		

\* Andere Seite: Beinstütze Standard

### Lower leg length in cm (without seat cushion)

Price

Lower leg length (LL) in cm  (available setting range 16 – 55 cm)



Foot support bar	Lower leg length [cm]					
	Selection foot bar by 0 cm front frame offset				Selection foot bar by 2 cm front frame offset	
	Leg support single panel, leg support individual		Elevating leg support left or right	Elevating leg support, pair	Leg support single panel, leg support individual	
	Foot plate Aluminium	Foot plate plastic			Foot plate Aluminium	Foot plate plastic
Short	33 - 43	38 - 43	30 - 39	28 - 39	32 - 40	34 - 42
Medium	37 - 48	43 - 48	–	–	37 - 45	–
Long	42 - 52	47 - 52	35 - 51	35 - 51	42 - 50	46 - 53
Extra long	47 - 55	–	–	–	47 - 54	–



# Active wheelchair *Avantgarde DV.*

## Master order form

### Choice of foot support bar (without seat cushion)

		Price
<input type="radio"/> MB40	Short	–
<input type="radio"/> MB41	Medium	–
<input type="radio"/> MB42	Long	–
<input type="radio"/> MB41	Extra long	–

### Leg support accessories

		Price
<input type="checkbox"/> MZ51	Lap belt with metal buckle	Add. charge
<input type="checkbox"/> MZ59	Pocket for mobile phone	Add. charge

### Spoke protectors

		Price
<input type="radio"/> MG90	Transparent	Add. charge
<input type="radio"/> MG100	Graphic decor	Add. charge
<input type="radio"/> MG101	Miracle decor	Add. charge
<input type="radio"/> MG102	Network decor	Add. charge
<input type="radio"/> MG103	Robots decor	Add. charge
<input type="radio"/> MG104	Spectrum decor	Add. charge

### Frame colours

		Price
<input type="radio"/> 001	Brilliant silver	–
<input type="radio"/> 002	Anthracite metallic	–
<input type="radio"/> 003	Jet black RAL 9005	–
<input type="radio"/> 004	Ice blue	–
<input type="radio"/> 005	Sparkle light orange	–
<input type="radio"/> 009	Vinho sparkle	–
<input type="radio"/> 010	Marine blue	–
<input type="radio"/> 012	Sky blue	–
<input type="radio"/> 014	Candy red	–
<input type="radio"/> 021	Sparkle granny-smith	–
<input type="radio"/> 023	White RAL 9016	–
<input type="radio"/> 027	Shimano matte	–
<input type="radio"/> 032	Pure metal blue	Add. charge
<input type="radio"/> 033	Pure metal green	Add. charge
<input type="radio"/> 034	Pure metal bronze	Add. charge
<input type="radio"/> 035	Pure metal red	Add. charge
<input type="radio"/> 048	Jet black matte (structure)	Add. charge
<input type="radio"/> SBK24817.000	Special colour to customer requirements	Add. charge

### Special accessories

		Price
<input type="checkbox"/> MZ52	Tool kit	Add. charge
<input type="checkbox"/> MZ58	Bag for drive wheels (for 2 drive wheels up to 28")	Add. charge
<input type="checkbox"/> MZ65	Active lighting	Add. charge
<input type="checkbox"/> MZ89	Protection for quick-release-axle	Add. charge
<input type="checkbox"/> MZ92	Wheelblades (Mini-skis for caster wheels)	Add. charge
<input type="checkbox"/> MZ93	Backpack	Add. charge



# Active wheelchair *Avantgarde DV.*

## Master order form

### Information

Caster wheels	3"	4"	5"	5,5"	6"
38,0 cm	●				
39,0 cm	●				
40,0 cm	●				
41,0 cm	●	●			
42,0 cm	●	●			
43,0 cm	●	●	●		
44,0 cm	●	●	●		
45,0 cm	●	●	●	●	●
46,0 cm	●	●	●	●	●
47,0 cm	●	●	●	●	●
48,0 cm		●	●	●	●
49,0 cm		●	●	●	●
50,0 cm		●	●	●	●
51,0 cm			●	●	●
52,0 cm			●	●	●
53,0 cm				●	●
54,0 cm					●

● possible combination

The values indicated are rounded values which have been theoretically determined. Max. deviation of 1 cm.

The anterior seat height depends on the selected wheel size, caster fork and mounting position. Please refer to the seat height table!

Specified without seat cushion at 0° seat inclination.

Horizontal position of the drive wheel attachment device: see sketch.


The posterior seat height (PSH) depends on the selected wheel size and mounting position in the drive wheel attachment device.

Please refer to the seat height table!

A manufacturing tolerance of 10 mm applies to all measurements and of 10% to specified weights.

### Wheelchair measurements

- 1 Anterior seat height (ASH) = ground to top of seat upholstery
- 2 Lower leg length (LL) = top of seat upholstery to top of foot plate (without seat cushion)
- 3 Posterior seat height (PSH) = ground to top of seat upholstery
- 4 Back support height (BH) = top of seat tube (cross-brace) to top of back support upholstery
- 5 Seat depth (SD) = front of seat upholstery to back support tube
- 6 Seat width (SW) = outside edge of frame tube to outside edge of frame tube





# Active wheelchair *Avantgarde DV.*

## Master order form

Horizontal drive wheel position A															
Seat depth [cm]	Drive wheel size 22"					Drive wheel size 24"					Drive wheel size 25"				
	Caster wheels														
	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"
36,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
38,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
40,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
42,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
44,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
46,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
48,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
50,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
52,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
54,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●

Horizontal drive wheel position B															
Seat depth [cm]	Drive wheel size 22"					Drive wheel size 24"					Drive wheel size 25"				
	Caster wheels														
	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"
36,0	●	●	●	●	●	●	●				●	●			
38,0	●	●	●	●	●	●	●				●	●			
40,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
42,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
44,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
46,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
48,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
50,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
52,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
54,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●

Horizontal drive wheel position C and D															
Seat depth [cm]	Drive wheel size 22"					Drive wheel size 24"					Drive wheel size 25"				
	Caster wheels														
	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"
36,0															
38,0															
40,0	●	●	●	●	●										
42,0	●	●	●	●	●	●	●	●	●	●					
44,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
46,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
48,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
50,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
52,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
54,0	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●



# Active wheelchair *Avantgarde DV.*

## Master order form

Horizontal drive wheel position E															
Seat depth [cm]	Drive wheel size 22"					Drive wheel size 24"					Drive wheel size 25"				
		Caster wheels													
	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"
36,0															
38,0															
40,0	●	●	●	●	●										
42,0	●	●	●	●	●	●	●	●							
44,0	●	●	●	●	●	●	●	●			●	●	●		
46,0	●	●	●	●	●	●	●	●			●	●	●		
48,0	●	●	●	●	●	●	●	●			●	●	●		
50,0	●	●	●	●	●	●	●	●			●	●	●		
52,0	●	●	●	●	●	●	●	●			●	●	●		
54,0	●	●	●	●	●	●	●	●			●	●	●		

Horizontal drive wheel position F and G																
Seat depth [cm]	Drive wheel size 22"					Drive wheel size 24"					Drive wheel size 25"					
		Caster wheels														
	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"	3"	4"	5"	5,5"	6"	
36,0																
38,0																
40,0																
42,0	●	●	●													
44,0	●	●	●			●	●	●								
46,0	●	●	●			●	●	●			●	●	●			
48,0	●	●	●			●	●	●			●	●	●			
50,0	●	●	●			●	●	●			●	●	●			
52,0	●	●	●			●	●	●			●	●	●			
54,0	●	●	●			●	●	●			●	●	●			

● possible position



# Active wheelchair *Avantgarde DV*.

## Master order form

Posterior seat height with wheelbase extension			
Axle position	Drive wheel size		
	22"	24"	25"
1	37 cm	39 cm	40 cm
2	40 cm	42 cm	42 cm
3	42 cm	44 cm	45 cm
4	46 cm	46 cm	48 cm
5	47 cm	49 cm	50 cm
Posterior seat height with shock absorber system			
Axle position	Drive wheel size		
	22"	24"	25"
1	40 cm	42 cm	43 cm
2	42 cm	44 cm	45 cm
3	44 cm	46 cm	47 cm
4	46 cm	48 cm	49 cm
5	48 cm	50 cm	51 cm
Posterior seat height			
Free mounting positions	Drive wheel size		
	22"	24"	25"
1	36 cm	38 cm	40 cm
2	37 cm	39 cm	41 cm
3	38 cm	40 cm	42 cm
4	39 cm	41 cm	43 cm
5	40 cm	42 cm	44 cm
6	41 cm	43 cm	45 cm
7	42 cm	44 cm	46 cm
8	43 cm	45 cm	47 cm
9	44 cm	46 cm	48 cm
10	45 cm	47 cm	49 cm
11	46 cm	48 cm	50 cm
12	47 cm	49 cm	51 cm

This table is for information only. The locking position may deviate in certain constellations. The locking position is selected such that the applicable rear seat height in the ordered constellation is reached as ordered.



# Active wheelchair *Avantgarde DV.*

## Master order form



**MA 55/56 Anti-tipper**



**MA 51/52 Tip-assist**




**MA 57/58 Crutch holder**



**MG83 One-arm drive with double handrim**



**MG10 Spinergy SPOX**




**MG06 Hollow rim**



**MG11 Spinergy LX**




**MG15 Infinity Ultralight**




**MG31 Right Run 1"**




**MG371 Marathon Plus Evolution 1.1"**



**MG28 Tetra release**



**MA 114 Drive wheel adapter, permanently welded**



**MG54 Handrim rubber-coated**



**MG65 Gekko**



**MG63 Curve L (eloxiert)**



**MG64 Curve L Grip**



## Active wheelchair *Avantgarde DV.*


### Master order form




**MF35/36/37 soft caster**  
4", 5", 6"



**MF30/31/32 Solid rubber,**  
black 4", 5", 6"



**MF04 Soft rubber 5,5"**



**MF05 Pneumatic tyre 6"**



**MD 56 Fold-down push**  
handles



**MD 55 Push handles, height-**  
adjustable/removable



**MH05 Scissor wheel lock,**  
aluminium



**MH 10 Wheel lock lever**  
extension



**MC 15 Seat upholstery**  
with storage compartments




**MB27 Leg support, elevating**



**MG80 Transparent**



**MG100 Graphics**



**MG101 Miracle**



**MG102 Network**



**MG103 Robots**



**MG104 Spectrum**