



When comfort comes first

Combining comfort, ease of handling and safety, the *Invacare SB 755* is perfectly suited to all home and long-term care environments.

The **SB 755** offers optimum comfort with an ideal seating position to minimise shear and friction. Carers can easily adjust the bed to meet individual client requirements.

Thanks to a unique technical feature, the **SB 755** can be adapted to suit taller clients by extending both the head and the foot end.

All these features provide real benefits to clients and carers in terms of comfort and general well-being.





Ease of handling

With the option to choose a bed that's specially adapted to withstand cleaning in a washing machine, the *Invacare SB 755* is now even easier to maintain.

Mounting and dismantling the **SB 755** has been made particularly easy to allow for greater ease of transport and storage. The bed can be stored on its base and rolled directly from van to bedroom. Assembly and transportation can be carried out quickly and with minimum effort by one person. No part of the bed weighs more than 25 kg and dismantlable scissors are available as an option.

Safety is key

As always at Invacare, safety is our key concern. The **SB 755** bed itself and side rails in both the standard bed length and the extendable version meet the strength demands of the **IEC 60601-2-52** safety standard for beds, without compromising the ergonomics or working conditions of the carer. Clients can count on this bed for total reliability and comfort – day or night.

Comfort is our concern

The **SB 755** comes with several features to ensure an excellent seating position and to reduce shear and friction during movement. These features include an auto-regression system, an optimised length of mattress platform sections and a leg

section that can be lowered below horizontal.


Invacare strives to cater for all clients and is delighted to offer an integrated length extension of 100 mm at the head and foot of this model to meet the needs of taller users.

Features and Options



Auto-regression function

Client comfort when activating backrest, minimising shear and friction.



Optimised dimensions of mattress platform according to anthropometrical measurement, for increased comfort.



Extendable mattress platform

Additional 100 mm at head and foot. Bed easily extends with 50, 100, 150 or 200 mm to suit client's needs.

Features and Options



NEW

Olivia bed end

Can be used in combination with wooden front panel. Available in light oak and beech.



Victoria bed end

Detachable panel. Available in beech, cherry, light oak and white beech.



Sophie bed end

Available in beech, cherry, light oak and white beech.



Emma bed end

Detachable - Available in beech, cherry, light oak and white beech.




Susanne bed end

Available in beech, cherry, white beech and light oak.



Anita bed end

Available in beech, cherry and light oak.



Britt V full-length side rail

Now available in 200, 210 and 220 cm length. Available in beech, white beech, cherry and light oak.



Line/ Extendable Line

Full-length aluminum side rails. Available in beech, white beech, cherry and light oak.



Scala 2 range

Collapsible side rail range in 3 different heights.




Wooden side panel

Available in beech, white beech, cherry and light oak.



Wooden side panel

Can be used in combination with Scala Basic Plus 2 side rail. Available in beech and light oak.



Mattress support

Mattress support with choice of aluminium or wooden slats.



Several handcontrols:

Invacare **Soft Control** or standard with "green" button function.



Lifting pole

The handle can be adjusted both in height and depth.



Swivel support handle

The swivel support handle provides good support for the client getting to a standing position.



Central brake

Available with a wheel diameter in 100, 125 or 150 mm.

Features and Options



Easy to store and transport

Bed parts, mattress and accessories can be stored directly on the base using special hooks, which means that storage and transport are simplified. The *Invacare SB 755* bed requires minimal storage space.

Technical data

SB 755							
	Width outside: 950 / 1000 mm Width inside: 850 / 900 mm	Length outside: 2130 mm Length inside: 2080 mm Mattress: 2000 mm	Height adjustable: 350-820 mm Total height: 1600-2070 mm	Mattress support dimensions: 800 - 260 - 260 - 680 mm	Under bed clearance: 170 mm	0° - 73° Backrest angle	0° - 28° Thigh angle
SB 755 Wide 1050 mm	Width outside: 1150 mm Width inside: 1050 mm	Length outside: 2130 mm Length inside: 2080 mm Mattress: 2000 mm	Height adjustable: 350 - 820 mm Total height: 1600 - 2070 mm	Mattress support dimensions: 800 - 260 - 260 - 680 mm	Under bed clearance: 170 mm	0° - 73° Backrest angle	0° - 28° Thigh angle
SB 755 Wide 1200 mm	Width outside: 1300 mm Width inside: 1200 mm	Length outside: 2130 mm Length inside: 2080 mm Mattress: 2000 mm	Height adjustable: 350 - 820 mm Total height: 1600 - 2070 mm	Mattress support dimensions: 800 - 260 - 260 - 680 mm	Under bed clearance: 170 mm	0° - 73° Backrest angle	0° - 28° Thigh angle


0° - 16° Minimum angle below horizontal	0° - 18° Full legrest lift (Angle below frame) 14° (Angle above frame)	25°	SWL 235 kg Max. user weight 200 kg	Heaviest part of the product: 25 kg	Washable in a washing machine	Bariatric use

Complies with IEC 60601-2-52

NB/ The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

Frame colour and wood work



Aluminium grey

Beech

Light oak

White beech

Cherry

© 2015 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

SB 755 - EU - 08/2015



Yes, you can.®